





SNOMED CT, NPU and LOINC integration for semantic interoperability of laboratory test results

Daniel Karlsson¹, Jay Kola², James Campbell³

- 1. National Board of Health and Welfare, Sweden
- 2. NHS Digital, UK
- 3. University of Nebraska Medical Center, US

Interoperability through SNOMED CT

- There are today a few laboratory medicine terminologies, e.g.
 - LOINC is used in US, AU, MY, and many more countries
 - NPU Terminology is used in Scandinavia, CZ
 - SNOMED CT is the lab terminology in UK
- Neither LOINC nor NPU are based on any ontologies



Interoperability through SNOMED CT

• Proposal to apply SNOMED CT to bridge between the three efforts



LOINC-SNOMED CT linking

- SNOMED International has since 2014 provided SNOMED CT concepts associated with LOINC parts and SNOMED CT expressions associated with LOINC terms
- Approx. 20 000 LOINC terms associated
- Project currently on hold



NPU-SNOMED CT pilot

- IUPAC project to determine feasibility of associating NPU Terminology to SNOMED CT
- Set of ~800 commonly used NPU terms
- Missing about 1/3 of elements/parts \rightarrow 3/4 of terms had 1 or more missing elements



SNOMED CT in OWL

- SNOMED International provides tools to render a SNOMED CT release in OWL format
 - <u>https://github.com/IHTSDO/snomed-owl-toolkit</u>
- LOINC expression association reference set can also be rendered as an OWL file



Lab terminologies on SNOMED CT

- Association between lab terminologies and SNOMED CT can facilitate
 - Querying
 - Use and re-use of lab data
 - Interoperability between users of different lab terminologies



)uonu	C	200	ovnroe	eion	1

Query (class expression)

'Observable entity (observable entity)' and

('Role group (attribute)' some ('Property (attribute)' some 'Quantity concentration (property) (qualifier value)')) and

('Role group (attribute)' **some** ('Inheres in (attribute)' **some** 'Urine (substance)')) **and**

('Role group (attribute)' some ('Component (attribute)' some 'Morphinan opioid (substance)'))

Levorphanol [Presence] in Urine (LOINC:3711-9)	0
Levorphanol [Presence] in Urine by Confirmatory method (LOINC:16213-1)	?
Levorphanol [Presence] in Urine by Screen method (LOINC:19516-4)	0
Morphine [Mass/volume] in Urine (LOINC:3831-5)	?
Morphine [Mass/volume] in Urine by Confirmatory method (LOINC:16251-1)	?
Morphine [Moles/volume] in Urine by Confirmatory method (LOINC:70210-0)	?
Morphine [Presence] in Urine (LOINC:3830-7)	?
Morphine [Presence] in Urine by Confirmatory method (LOINC:16196-8)	?
Morphine [Presence] in Urine by Screen method (LOINC:19597-4)	?
• NPU04536 U—Nalbuphine; arb.c.(proc.) = ?	?
•NPU04591 U—Oxycodone; arb.c.(proc.) = ?	0
•NPU08990 U—Morphine+analogue; arb.c.(proc.) = ?	0
^e NPU23881 U—Ethylmorphine; mass c. = ? μg/L	0
[●] NPU24866 U—Morphine; mass c. = ? μg/L	0
•NPU26819 U—Morphine; arb.c.(confirm; proc.) = ?	0
•NPU27984 U—Morphine+analogue; arb.c.(confirm; proc.) = ?	0
^e NPU28000 U—Oxycodone; mass c. = ? μg/L	?
Nalbuphine [Mass/volume] in Urine (LOINC:3840-6)	0
Nalbuphine [Mass/volume] in Urine by Confirmatory method (LOINC:20551-8)	0
Nalbuphine [Presence] in Urine (LOINC:3839-8)	?
Nalbuphine [Presence] in Urine by Confirmatory method (LOINC:16218-0)	?
Nalbuphine [Presence] in Urine by Screen method (LOINC:19605-5)	0
Oxycodone [Mass/volume] in Urine (LOINC:11246-6)	0
Oxycodone [Mass/volume] in Urine by Confirmatory method (LOINC:16249-5)	0
Oxycodone [Moles/volume] in Urine by Confirmatory method (LOINC:70215-9)	?
Oxycodone [Presence] in Urine (LOINC:10998-3)	2

Sodium in (Serum or) Plasma

Description: 'Sodium [Moles/volume] in Serum or Plasma (LOINC:2951-2)'

Equivalent To 🕀

```
Clinical chemistry observable (LOINC grouper)'
and ('Role group (attribute)' some ('Component (attribute)' some 'Sodium electrolyte (substance)'))
and ('Role group (attribute)' some ('Property (attribute)' some 'Substance concentration (property) (qualifier val
and ('Role group (attribute)' some ('Scale type (attribute)' some 'Quantitative (qualifier value)'))
and ('Role group (attribute)' some ('Time aspect (attribute)' some 'Single point in time (qualifier value)'))
and ('Role group (attribute)' some ('Inheres in (attribute)' some 'Plasma (substance)'))
and ('Role group (attribute)' some ('Direct site (attribute)' some 'Acellular blood (serum or plasma) specimen
(specimen)'))
```

SubClass Of 🕀

```
'Serum or plasma clinical chemistry (observable entity)'
```

```
Sodium [Moles/volume] in Body fluid (LOINC:2950-4)
```


BINDH02420 D. Codium ion: cubet a - 2 mmol/L!

Description: NPU03429 P—Sodium ion; subst.c. = ? mmol/L

Equivalent To 💽

Observable entity (observable entity)' and ('Role group (attribute)' some ('Component (attribute)' some 'Sodium electrolyte (substance)')) and ('Role group (attribute)' some ('Property (attribute)' some 'Substance concentration (property) (qualifier value)')) and ('Role group (attribute)' some ('Inheres in (attribute)' some 'Plasma (substance)'))





('Role group (attribute)' some ('Property (attribute)' some 'Susceptibility (property) (qualifier value)')) and

('Role group (attribute)' some ('Inheres in (attribute)' some 'Bacterium (organism)'))

SubClassOf

Role group (attribute)' some ('Property (attribute)' some 'Threshold mass concentration (property) (qualifier value)')

Socialstyrelsen

Conclusions

 SNOMED CT and the Observables model can be used as an interlingua for lab terminology interoperability



Thanks!

